

Loupes & Lights




enova®

WE PRIDE OURSELVES ON...



ERGONOMICS

Creating all day
comfort



DETAIL

Quality to the
highest standard



RESOLUTION

Highest resolution
and the cleanest
viewing experience



CUSTOMER SERVICE


Five star service
guaranteed



INTENSITY

World's brightest LED
surgical headlight

Enova Vizix Loupes



Modern Black



Silver Snow



Purple Rain



Lake Blue



Piper Green

**MADE IN
USA** 



2.5x



3.0x



3.5x

Vizix Frame Options



Flak 2.0 XL



Latch-Alpha



MainLink



Curved



Titanium

Axis Light Options



Activ Battery Option

- 5+ hour run-time
- Comes with 2 batteries
- Capacitive touch technology
- 3 levels of intensity
- Silicone Anti-slip Neck Band
- Color Customizable



Standard PP1.5 Battery Option

- 10+ hour run-time
- 12 levels of intensity
- Attaches to the hip

Color Customizable
& Optional Yellow Filter

Clips on magnetically

Aircraft Aluminum LED Housing

24,000 Lux



Modern Black



Silver Snow



Purple Rain



Lake Blue



Piper Green

Atom Light Options



- 25,000 Lux
- Wireless
- Clips on magnetically
- Three rechargeable batteries
- 2 hour run-time
- Capacitive Touch Technology
- 2 Levels of Intensity
- Color Temperature 5700 Kelvin
- Weight: .25 oz / 7g



MINNESOTA MADE

● 1717 Broadway Street NE ● Minneapolis, MN 55413 ● 651.236.8857 ● goenova.com ●