



Quasar Cool (XLT-150A) Coaxial LED Surgical Headlight

Includes:



charging station



carrying case



2 PP2 batteries



battery cables



battery belt

Specifications:

Light Source: LED refractive technology

Lamp Illumination: **150,000 - 300,000 LUX**

Color Temperature: **6100K**

Illumination Spot Size: adjustable from 2" to 5" (5.08cm to 12.7cm)

Headlight Construction: all aluminum

Headlight Weight: 6oz (170g)

Headband Construction: rigid plastic, double removable foam pads

LED Warranty: 5 years

Battery: rechargeable lithium-ion (charger included)

Battery Life: 10 to 24 hours with the PP4 (4-cell) battery;
5 to 10 hours with PP2 (2-cell battery);
4 to 8 hours with Freestyle Battery

Battery Pack Dimensions: PP4 is 3 X 4 X 1"

Battery Pack Weight: PP4 is 9 oz/255gm

Battery Warranty: 1 year

Free evaluation goenova.com/evaluation-program/