

# Quasar Cool (XLT-150A)




## Enova's Most Versatile LED Surgical Headlight

150,000 - 300,000 LUX

Adjustable spot 2" - 5" / 5cm - 13cm

Recommended for physicians  
performing procedures for:

- Plastic, Cosmetic, & Reconstructive
- Oral and Maxillofacial
- Dental – General and Endo
- General Surgery
- EM/ER
- ENT
- OB/GYN
- Neuro
- Vascular
- Vet – General & Small Animal



### REFRACTIVE LED

Brightness (LUX)

150,000 LUX

Color Temperature (Kelvin)

6100 K

Color Rendering Index (out of 100)

CRI 75

### RECHARGABLE LITHIUM-ION

Freestyle (2-cell)

4 - 8 Hours

PP2 (2-cell) Battery

5 - 10 Hours

PP4 (4-cell) Battery

9 - 24 Hours