

QUANTUM™ COOL LED HEADLIGHT



Adjustable light angle
for ergonomic wear

Adjust fit from over
and behind the
head for comfort
and no-slip wear

Durable Aluminum Housing

All comfort pads
treated with Aegis®

Adjustable spot

THE WORLD'S BRIGHTEST LED HEADLIGHT

The illumination spot of the **Quantum Cool LED Headlight** features clear, cool, white light that is uniform and consistent. This headlight offers the adjustability that makes it our most versatile and popular surgical headlight.

250,000 - 500,000 LUX

Utilizing three LEDs under one bezel, the Quantum Cool headlight packs the brightest punch on the market, offering 250,000-500,000 lux.*

COMFORTABLE AND ERGONOMIC

Weighing only 12oz, the Quantum headlight is one of the lightest, most comfortable, and durable surgical headlights.

INFECTION PREVENTION

Enova headlights do not require fans, all headbands are treated with antimicrobial coating, and easy-to-clean light housings simplify infection prevention and control.

*Lux is variable based on preferred working distance.

TECHNICAL DATA

Light Source: LED refractive technology

Lamp Illumination: 250,000 (at 14") - 500,000 (at 10") lux

Color Temperature: 6,100K

Illumination Spot Size: Adjustable from 2" to 5" (5.08cm to 12.7cm)

Headlight Construction: All Aluminum

Headlight Weight: 12oz (340g)

Headband Construction: Nylon Copolymer, Replaceable Comfort Pads

LED Warranty: 5 years

Battery: Rechargeable Lithium Ion (base station charger included)

Battery Warranty: 1 year

Battery Pack Options	PP4
Battery Run-Time	10-24 hours
Battery Pack Dimensions	3"x4"x1"
Battery Pack Weight	9 oz (255g)

DESIGNED TO MEET YOUR NEEDS

Enova offers a range of the brightest, lightest surgical headlamps so you can select the exact specifications to fit your needs. Try an Enova headlight with our free evaluation program! Visit www.goenova.com to learn more.



To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road STE 100
Minneapolis, MN 55427

