

# VIZIX<sup>®</sup> LOUPES

## Prismatic



Every Vizix frame comes with a magnetic adapter to securely house clip-on lights or cameras.

Built with all aluminum housing for easy cleaning and greater durability than plastic loupes.

All-glass telescopes for high quality optical resolution, clarity, and sharpness.

Prescription-friendly telescopes and optional carrier lens to create a precision fit unique to your eyes.

## High Quality Magnification, Made for You

Experience crisp, clear vision with Enova's **Vizix<sup>®</sup> Prismatic Loupe**. With magnifications ranging from 3.5x to 5.5x, you can select the perfect level of magnification for your specific needs - ensuring maximum performance and comfort during any professional task throughout the day! Experience pristine clarity - get ready for crystal-clear viewing now!

### MAGNIFICATION

#### 3.5x - 5.5x

Select the level of magnification that best fits your needs. Our experienced representatives can help you identify a magnification level that is right for you.

### CUSTOM

#### MEASUREMENTS

Elevate the way you work with a custom-made pair of loupes that suit your exact needs - whether it's working distance, declination angle or lens prescription. Enhance performance and be on top of your game!

### FUNCTIONAL

#### FRAMES

Enova frames provide the complete package: modern style and superior durability. Comfort, ergonomics, and practicality are all seamlessly integrated into an attractive design.

# TECHNICAL DATA

**Ocular Construction:** Aluminum housing and all-glass telescope

**Ocular Warranty:** Lifetime limited warranty on manufacturer's defects

## \*Ocular Weight

Ocular Model	3.5x	4.5x	5.5x
Prismatic P-Series	44.6g	50.4g	55.8g

\*Weight of pair with zero power, glass diopters.

## Frame Pairing Specifications

Model	*Frame Weight (g)	Material	Lens Height (mm)	Frame Width (mm)	Bridge Width (mm)	Temple Length (mm)
Airon Radius	15.8	Aluminum	42.8	149.2	11.7	130.5
Airon Standard	17.9	Aluminum	43.0	149.3	12.2	142.0
Airon Small	15.1	Aluminum	40.7	137.7	8.8	136.0
Airon Sport	13.8	Aluminum	50.2	145.7	8.3	135.0

\*Frame weighed without carriers or oculars

## DESIGNED TO MEET YOUR NEEDS

Take your loupes to the next level with Enova's selection of accessories! Our 4K cameras, clip-on lights, and more are sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success. Visit [www.goenova.com](http://www.goenova.com) or email customer service at [sales@goenova.com](mailto:sales@goenova.com) to learn more.



### Explore More Loupes



To learn more about how Enova can help you see better, so you can perform better, visit [www.goenova.com](http://www.goenova.com), call 651.236.8857, or email [sales@goenova.com](mailto:sales@goenova.com).

### World Headquarters

Enova Illumination  
6305 Sandburg Road STE 100  
Minneapolis, MN 55427

M22008A

© 2023 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

