

# VIZIX® LOUPES

## Angled Prismatic



Every Vizix frame comes with a magnetic adapter to securely house clip-on lights or cameras.

Provides unparalleled range of visibility through the lenses

Built with all aluminum housing for easy cleaning and greater durability than plastic loupes.

All-glass telescopes for high quality optical resolution, clarity, and sharpness.

Prescription-friendly telescopes and optional carrier lens to create a precision fit unique to your eyes.

## Lightweight, Ergonomic, Built for You

The Enova **Vizix® Angled Prismatic Loupe** is the ultimate accessory to enhance your practice. Modern engineering, found in Enova Vizix Loupes, places the telescope lower on the carrier lens for a greater range of vision. When used properly, in the right setting, it will help decrease fatigue and pain associated with repetitive cervical spine movement. The Angled Prism is available in 3.5x and 4.5x coming soon! All carrier frames have been hand-selected by experts for ergonomics, comfort, and durability. Select from Enova's Airon® family of frames for a stylish presentation.

### GREATER RANGE OF VISION

Strategically placed telescopes provide an expansive range of vision, allowing you to easily work around equipment without having to take off your loupes.

### LIGHTER & MORE COMPACT

The telescope is designed to be lightweight, without sacrificing technology that supports decreasing stress on your neck, shoulders, and lower back for long-term body health.

### FUNCTIONAL FRAMES

Enova frames provide the complete package: modern style and superior durability. Comfort, ergonomics, and practicality are all seamlessly integrated into an attractive design.

# TECHNICAL DATA

**Ocular Construction:** Aluminum housing and all-glass telescope

**Ocular Warranty:** Lifetime limited warranty on manufacturer's defects

## \*Ocular Weight

Magnification	3.5x
Angled Prismatic	32.8g

*\*Weight of pair with zero power, glass diopters.*

## Frame Pairing Specifications

Model	*Frame Weight (g)	Material	Lens Height (mm)	Frame Width (mm)	Bridge Width (mm)	Temple Length (mm)
Airon Radius	15.8	Aluminum	42.8	149.2	11.7	130.5
Airon Standard	17.9	Aluminum	43.0	149.3	12.2	142.0
Airon Small	15.1	Aluminum	40.7	137.7	8.8	136.0
Airon Sport	13.8	Aluminum	50.2	145.7	8.3	135.0

*\*Frame weighed without carriers or oculars*

## DESIGNED TO MEET YOUR NEEDS

Take your loupes to the next level with Enova's selection of accessories! Our 4K cameras, clip-on lights, and more are sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success. Visit [www.goenova.com](http://www.goenova.com) or email customer service at [sales@goenova.com](mailto:sales@goenova.com) to learn more.



**4K Auto-focus Cameras**  
(proCam XS)



**Axis™ LED Light System**  
30,000 lux



**Atom Plus Wireless LED Light System**  
26,000 lux

### Explore More Loupes



**Prism Loupes**



**Galilean Loupes**

To learn more about how Enova can help you see better, so you can perform better, visit [www.goenova.com](http://www.goenova.com), call 651.236.8857, or email [sales@goenova.com](mailto:sales@goenova.com).

### World Headquarters

Enova Illumination  
6305 Sandburg Road STE 100  
Minneapolis, MN 55427

