

 enova®



**SEE BETTER
PERFORM
BETTER**
ZUMAX® Microscopes

zumax

REVOLUTIONIZE YOUR PRACTICE

WITH A ZUMAX® MICROSCOPE

Magnification using the latest in microscope technology gives professionals like you the opportunity to offer a higher standard of care in dentistry. Enova Illumination is proud to be the exclusive U.S. distributor of Zumax® microscopes, bringing you the highest quality and latest innovations available in the USA.



ENHANCED MAGNIFICATION

Increased precision in all aspects of patient care - from diagnosis and treatment, to efficiency and patient outcomes.



INTEGRATED 4K CAMERA

Connect to a monitor or your network for the ability to stream and save high definition photos and videos. Edit videos immediately for training and sharing with consulting physicians. Display real-time imaging for increased patient education.



SUPERIOR ERGONOMICS

Reduce neck and back injuries that can limit a career by not relying on loupes and insufficient illumination. Improve your posture by utilizing heads-up magnification with a monitor. 3D visualization is optional.



DOCUMENTATION AND TRAINING

Microscopes allow you to easily show, discuss and document treatments and procedures, both with patients and other specialists.



HIGH QUALITY OPTICS

On all Zumax® microscopes, all lenses are multi-coated, anti-reflective, apochromatic, and made from hand-selected German glass.



RETURN ON INVESTMENT

Elevate your practice with Zumax® Microscopes, empowering you to confidently perform complex procedures with enhanced visual precision!

Zumax® Dental Microscopes

3200 R2/PRO

FEATURES

- **Joystick-Controlled Continuous Zoom:**
Magnification ranging from ~1-20x
AutoFocus is standard
- **Electromagnetic Clutch:**
One-click magnetic clutch for ease of
positioning and smooth movement
- **Integrated 4K Camera:** Seamlessly save,
record and stream images/videos to your
network and monitor(s)
*3D 4K camera optional
- **Integrated LED Illumination:**
≥105,000 lux at 250mm working distance,
multiple spot sizes and light and laser
filters included. LED lifetime: 60,000 hours
- **Integrated VarioScope:**
Working range 200mm to 450mm



**Motorized Continuous
Zoom System**



**Electromagnetic
Clutch**



**Integrated 4K
Camera**

Zumax® Dental Microscopes

2380 R2

FEATURES

- **Continuous Zoom Magnification:** Ranging from ~1-20x
- **FreeFlow balancing arm:** Free flowing, ultra smooth movement to simplify positioning
- **Integrated 4K Camera:** Seamlessly save, record and stream images/videos to your network and monitor(s)
*3D 4K camera optional
- **Integrated LED Illumination:**
≥100,000 lux at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **Integrated VarioScope:**
Variable working distance 200mm-450mm



**Continuous
Zoom System**



**VarioScope Fine
Focus**



**Integrated 4K
Camera**

Zumax® Dental Microscopes

2050

FEATURES

- **Continuous Zoom Magnification:**
Ranges from ~1-20x
- **ErgoPendulum Arm:**
Binoculars remain stationary as microscope pivots right and left
- **Integrated LED Illumination:**
≥75,000 LUX at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **4K Camera Module (Optional):**
Seamlessly save and display images and video
- **Tabletop Mount:**
Available option



**Continuous
Zoom System**



Integrated LED

Zumax® Dental Microscopes

2350

FEATURES

- **6-Step Magnification:**
Ranging from ~1-24x
- **Integrated LED Illumination:**
≥65,000 LUX at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **Integrated HD Camera (Optional):** Seamlessly connects to your network
- **Integrated VarioScope:**
Variable working distance 190mm-300mm



Integrated LED



VarioScope Fine Focus

ACCESSORIES

BROKEN INSTRUMENT REMOVAL KIT

Instrument separation is one of the most common complications in the root canal treatment process. The innovative design of the Zumax® Broken Instrument Removal Kit makes instrument removal quick and easy.

Kit Includes:

Extractor Handle and Tip



Extractor Tip Options



- 3 Micro-Trephine Burs
- 1 Extractor Handle
- 3 Extractors
- 1 Foreign Forceps

MICRO DENTAL TWEEZER

The Zumax® Micro Dental Tweezer is designed for accurate foreign material removal - no matter the angle or complexity.

Kit Includes:

- 1 360 Tweezer Head
- 1 Handle
- 1 Cleaning Adapter



MOBILE ADAPTERS

Connect mobile devices to Zumax® microscopes to allow you to instantly share your images through live-streaming, text messages, or email.



Contact your Enova sales rep today to learn more about Zumax microscopes and accessories.
Visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

ENOVA SATISFACTION GUARANTEE

Use our products and love them. If not, we'll take them back. It's that simple.
For full warranty specifications visit www.goenova.com/product-support

To learn more about how Enova can help you see better, so you can perform better,
visit www.goenova.com, call 651.236.8857, or email sales@goenova.com

World Headquarters

Enova Illumination
6305 Sandburg Road, STE 100
Minneapolis, MN 55427

ISO 13485: 2016 Certified QMS
MDR Class 1 Certified

PLACEHOLDER FOR DOCUMENT NUMBER

© 2023 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Oakley® and Flak® are trademarks of Oakley, Inc.
Zumax® is a registered trademark of Zumax Medical Co., LTD

