

AXIS™ series

DENTAL LIGHT



ADVANCED COMFORT AND VISUALIZATION:

Enova developed these **Axis™ LED Magnetic Lights** to be lightweight, an ergonomic illumination solution compatible with all Enova loupes. Either light can be paired with the Activ neck battery or PP1.5 hip battery to provide variable light intensity, and maximize magnification for the field of view. When used in conjunction with Enova's lightweight Vizix® loupes, clinicians benefit from superior illumination, visualization, and comfort, allowing them to perform with ease and confidence.

AXIS ELITE: LIGHT WEIGHT

The new Axis Elite is one of the lightest dental lights on the market. The pivot connection allows for personalized alignment with field of view.

AXIS STANDARD: PERSONALIZATION

Choose from five different distinct color caps to personalize your light. Available in silver, black, blue, green and purple. Activ™ Battery caps also come in black, green, blue and purple to complete the look.

OPTIONAL HD GLASS COMPOSITE FILTER

Removes the light wavelength that tends to speed up the curing of composite materials specific to the dental field. When not needed, simply flip the filter up for pure LED light.

*Lux is variable based on preferred working distance

TECHNICAL DATA

	Axis Elite	Axis Standard	Atom Plus
Image & Color Options			
Light Source	High Performance LED	High Performance LED	LED Refractive Technology
Brightness	40,000 Lux @ 10 in	30,000 Lux @ 10 in	26,500 Lux @ 10 in
Color Temp.	5,700K	6,100K	6,000K
Light Spot Size	3.25 in @ 12 in	4 in @ 12 in	4 in @ 12 in
Housing Material	Aluminum	Aluminum	Aluminum
Compatible Battery	Activ Elite (Neck Battery) PP1.5 Elite (Hip Battery) <small>Not interchangeable</small>	Activ (Neck Battery) PP1.5 (Hip Battery)	Wireless Magnetic Battery
Battery Run-Time	Activ Elite: 5-8 hours PP1.5 Elite: 10-20 hours	Activ: 5-8 hours PP1.5: 10-20 hours	3+ hours
Light Weight	0.28oz, 7.8g (w/o filter)	0.62oz, 17.6g (w/o filter)	0.44oz, 12.6g (w/o filter)
Battery Weight	Activ Elite: 5.3oz (150g) PP1.5 Elite: 6.5oz (185g)	Activ: 5.3oz (150g) PP1.5: 6.5oz (185g)	0.8oz (23g)
Battery Size	Activ Elite: 1×1 in, 21 in long PP1.5 Elite: 4×2.5×1 in	Activ: 1×1 in, 21 in long PP1.5: 4×2.5×1 in	2.2×1.1×0.6 in

	Axis Elite	Axis Standard	Atom Plus
FEATURES			
Compatible adapters to fit most manufacturers' loupes	●	●	●
Touch control - Brightness	●	●	●
(Not available with PP1.5 or PP1.5 Elite batteries)			
PACKAGE: DEFAULT vs OPTION			
UV filter (non-detachable)	●	●	●
No UV filter	●	●	●
Detachable light connector	●	—	—
Wired battery (x2 selected type(s), Rechargeable Li-Ion)	●	●	—
Wireless battery (x2, Rechargeable Li-Ion)	—	—	●
Single battery charger	●	●	—
Double battery charger	—	—	●
Custom engraving on the case	●	●	—
WARRANTY			
LED covering manufacturer defects - 3 years	●	●	●
Batteries and accessories - 1 year	●	●	●



Connect with your sales representative now to experience a whole new light.

Enova Illumination

O: 651.236.8857

E: sales@goenova.com

W: goenova.com

A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20218 Rev. C

Copyright © 2025 Enova Illumination, Inc. All rights reserved.
All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.