

# ATOM™ PLUS

## WIRELESS DENTAL LIGHT



## A TETHER-FREE LIGHT SOLUTION:

Experience the freedom of a wireless, tangle-free work environment with **Atom Plus Wireless LED Light**. This lightweight, versatile solution provides a bright, white spotlight that reaches 26,500 lux at maximum intensity. With magnetic clip batteries that attach easily and wirelessly without cables, you'll have a streamlined tangle-free working environment. The Atom Plus wireless light is ideal for professionals who need a powerful lighting option for clinical or minor surgical procedures. The package comes with two rechargeable batteries, a charger, detachable filter, and screwdriver- all packed in a carrying case.

### VERSATILE AND ADAPTABLE

Recommended for clinical professionals performing procedures for:

- Facial Plastics
- Oral and Maxillofacial
- General Dentistry
- Cosmetic Dentistry

### MAGNETIC LOUPE ADAPTERS

With magnetic loupe adapters, quickly and easily fit your existing loupes or most major brands with the perfect accessory. The reliable magnet attachment will provide all day hold.

### TOUCH ACTIVATED CONTROL

Touch activation powers the light and controls three levels of brightness intensity, providing you flexibility during procedures or exams.

\*Lux is variable based on preferred working distance

# TECHNICAL DATA

|                       | Axis Elite  | Axis Standard  | Atom Plus   |
|-----------------------|---|--|---|
| Image & Color Options |  |  |  |
| Light Source          | High Performance LED  | High Performance LED   | LED Refractive Technology   |
| Brightness            | 40,000 Lux @ 10 in  | 30,000 Lux @ 10 in   | 26,500 Lux @ 10 in  |
| Color Temp.           | 5,700K  | 6,100K   | 6,000K  |
| Light Spot Size       | 3.25 in @ 12 in   | 4 in @ 12 in   | 4 in @ 12 in  |
| Housing Material      | Aluminum  | Aluminum   | Aluminum  |
| Compatible Battery    | Activ Elite (Neck Battery)<br>PP1.5 Elite (Hip Battery)                           | Activ (Neck Battery)<br>PP1.5 (Hip Battery)  | Wireless Magnetic Battery   |
| Battery Run-Time      | Activ Elite: 5-8 hours<br>PP1.5 Elite: 10-20 hours                                | Activ: 5-8 hours<br>PP1.5: 10-20 hours   | 3+ hours  |
| Light Weight          | 0.28oz, 7.8g (w/o filter)   | 0.62oz, 17.6g (w/o filter)   | 0.44oz, 12.6g (w/o filter)  |
| Battery Weight        | Activ Elite: 5.3oz (150g)<br>PP1.5 Elite: 6.5oz (185g)                            | Activ: 5.3oz (150g)<br>PP1.5: 6.5oz (185g)   | 0.8oz (23g)   |
| Battery Size          | Activ Elite: 1x1 in, 21 in long<br>PP1.5 Elite: 4x2.5x1 in                        | Activ: 1x1 in, 21 in long<br>PP1.5: 4x2.5x1 in                                     | 2.2x1.1x0.6 in  |

|  | Axis Elite | Axis Standard | Atom Plus |
|--|------------|---------------|-----------|
| <b>FEATURES</b>  |            |               |           |
| Compatible adapters to fit most manufacturers' loupes    | ●          | ●             | ●         |
| Touch control - Brightness                               | ●          | ●             | ●         |
| (Not available with PP1.5 or PP1.5 Elite batteries)      |            |               |           |
| <b>PACKAGE: DEFAULT vs OPTION</b>                        |            |               |           |
| UV filter (non-detachable)                               | ●          | ●             | ●         |
| No UV filter   | ●          | ●             | ●         |
| Detachable light connector                               | ●          | —             | —         |
| Wired battery (x2 selected type(s), Rechargeable Li-Ion) | ●          | ●             | —         |
| Wireless battery (x2, Rechargeable Li-Ion)               | —          | —             | ●         |
| Single battery charger                                   | ●          | ●             | —         |
| Double battery charger                                   | —          | —             | ●         |
| Custom engraving on the case                             | ●          | ●             | —         |
| <b>WARRANTY</b>  |            |               |           |
| LED covering manufacturer defects - 3 years              | ●          | ●             | ●         |
| Batteries and accessories - 1 year                       | ●          | ●             | ●         |



Connect with your sales representative now to experience a whole new light.

## Enova Illumination

O: 651.236.8857

E: sales@goenova.com

W: goenova.com

A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20445 Rev. A

Copyright © 2025 Enova Illumination, Inc. All rights reserved.  
All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.