

PRISMATIC

VIZIX® LOUPES



Built with all aluminum housing for easy cleaning and greater durability than plastic loupes

HIGH QUALITY MAGNIFICATION:

Experience crisp, clear vision with Enova's **Vizix® Prismatic Loupe**. With magnifications ranging from 3.5x to 5.5x, you can select the perfect level of magnification for your specific needs - ensuring maximum performance and comfort during any professional task throughout the day! Experience pristine clarity - get ready for crystal-clear viewing now!

MAGNIFICATION 3.5x, 4.5x, 5.5x

Select the level of magnification that best fits your needs. Our experienced representatives can help you identify a magnification level that is right for you.

CUSTOM MEASUREMENTS







Elevate the way you work with custom-made loupes that suit your exact needs. Enhance your performance and be on top of your game!

FUNCTIONAL FRAMES

Enova frames provide the complete package: modern style and superior durability. Comfort, ergonomics, and practicality are all seamlessly integrated into an attractive design.











OCULAR SPECIFICATIONS





*Oculars weighted as a pair without frame

	Galilean	Prismatic	Ergo Elite
Image			
Ocular Color	 *Some colors may not be available in 3.5x		
Magnification	2.5x / 3.0x / 3.5x	3.5x / 4.5x / 5.5x	3.5x / 4.5x / 5.5x
Ocular Weight	2.5x: 28.8g, 3.0x: 34.4g, 3.5x: 36.4g	3.5x: 44.6g, 4.5x: 50.4g, 5.5x: 55.8g	3.5x: 32.8g, 4.5x: 35.7g, 5.5x: 35.2g
Housing Material	Aircraft-grade Aluminum	Aircraft-grade Aluminum	Aluminum, Plastic
Ocular Lens Material	Optical glass with anti-scratch & anti-reflective coating		

FRAME SPECIFICATIONS


*Frames weighted without lens or oculars

Enova Brand	Airon™ Standard	Airon™ Small	Airon™ Radius	Airon™ Tortoise	Airon™ Sport
Image					
Frame Color					
Frame Weight	17.9g	15.1g	15.8g	20.8g	16.7g
Frame Material	Aluminum	Aluminum	Aluminum	Acetate, Aluminum	Aluminum
Carrier Rx Material	High Index Trivex® with anti-reflective coating				Nylon

Oakley® Brand	Flak® 2.0 XL	CenterBoard
Image		
Frame Color		
Frame Weight	16.7g	12.9g
Frame Material	O-Matter, Unobtanium	Plastic
Carrier Rx Material	High Index Trivex® with anti-reflective coating	

EXPLORE MORE !


Take your loupes to the next level with Enova's selection of accessories! Our magnetic lights, 4K cameras, and more are sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success.



Axis™ Elite
LED Light System



Axis™ Standard
LED Light System



Atom™ Plus Wireless
LED Light System



To learn more about how Enova can help you see better, so you can perform better, visit our website or contact to us.

Enova Illumination

O: 651.236.8857

E: sales@goenova.com

W: www.goenova.com

A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20620 Rev. C

Copyright © 2025 Enova Illumination, Inc. All rights reserved.
All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.