

QUBIT™ PRO

LED HEADLIGHT



All Comfort Pads Treated with MICROBAN®

Adjustable Headband for Customizable Fit

Adjustable Links for Ergonomic Wear

Durable Aluminum Housing

Fixed Spot Size



Beam Adjustment Filter (optional)

Doubles the otherwise fixed spotlight size.

THE MOST ECONOMICAL LED HEADLIGHT

The illumination spot of the **Qubit Pro LED Headlight** features clear, cool, white light that is uniform and consistent. This headlight offers a lightweight and economical way to bring LED light into any practice. Attach the the Beam Adjustment Filter to make the headlight semi-adjustable.

55,000 - 120,000 LUX

Ranging from 55,000-120,000 lux* - depending on working distance - of pure white homogeneous light for extreme surgical field illumination.

COMFORTABLE AND ERGONOMIC

Designed for ultimate adjustability, comfort, and durability, the Quantum Cool headlight adjusts to your exact fit, desired working distance and light angles.

BEAM ADJUSTMENT

A simple flip filter allows the Qubit Pro to have a semi-adjustable light that doubles the spot size. For smaller, more focused light, simply flip the filter off the lens.

*Lux is variable based on preferred working distance

TECHNICAL DATA

Model Name	Quest	Quantum Cool	Quasar Cool	Qubit Warm	Qubit Pro
Light Source	High-powered LED	LED refractive technology			
Brightness @ 10"	240,000 Lux	500,000 Lux	300,000 Lux	160,000 Lux	230,000 Lux
Brightness @ 14"	125,000 Lux	250,000 Lux	150,000 Lux	80,000 Lux	115,000 Lux
Color Temp.	6,200K	6,500K	6,500K	4,500K	4,500K
Spot Size @ 14"	2.8-4.8 in	2-5 in	2-5 in	2 in	2-5 in
Headlight Housing	Thermoplastic, Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
Headlight Weight	10.7oz/303g	12.0oz/341g	6.5oz/183g	6.6oz/186g	6.6oz/188g
Compatible Battery	PP4+ (Plus)	PP4	PP2 or PP4		
Battery Run-Time	PP4+ (Plus): 4-18 hrs	PP4: 5-10 hrs	PP2: 5-10 hrs, PP4: 10-24 hrs		
Battery Weight	PP4+ (Plus): 10.3oz/292g	PP4: 10.6oz/300g	PP2: 6.5oz/184g, PP4: 10.6oz/300g		
Battery Size	PP4+ (Plus): 3.2×4×1 in	PP4: 3.2×4×1 in	PP2: 2.3×4×1 in, PP4: 3.2×4×1 in		

MODEL INFORMATION

	Quest	Quantum Cool	Quasar Cool	Qubit Warm	Qubi Pro
FEATURES					
Adjustable spot size	●	●	●	–	–
Pivot bracket	●	●	●	●	●
Hybrid thermal management using fan	●	–	–	–	–
Upgrade battery	–	–	●	●	●
Optional UV filter	–	–	●	●	●
PACKAGE: DEFAULT vs OPTION					
Single ratchet	●	●	●	●	●
Single ratchet w/occipital	●	●	●	●	●
Double ratchet	●	●	●	●	●
Double ratchet w/occipital	●	●	●	●	●
PP2 battery (Rechargeable Li-Ion)	–	–	●	●	●
PP4 battery (Rechargeable Li-Ion)	–	●	●	●	●
PP4 Plus battery (Rechargeable Li-Ion)	●	–	–	–	–
WARRANTY					
Headlamp covering manufacturer defects - 5 years	–	●	●	●	●
Headlamp covering manufacturer defects - 1 year	●	–	–	–	–
Headlamp pigtail - 1 year	●	●	●	●	●
Batteries and accessories - 1 year	●	●	●	●	●



Test the brilliance of the Enova headlight through our complimentary evaluation program.
Connect with your sales representative now to experience surgery in a whole new light.

Enova Illumination

O: 651.236.8857

E: sales@goenova.com

W: goenova.com

A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20623 Rev. B

Copyright © 2025 Enova Illumination, Inc. All rights reserved.

All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

Microban® is a registered trademark of Microban Products Company.

Antimicrobial Technology does not protect the user against disease-causing pathogens.