

# ATOM PLUS WIRELESS LED LIGHT



## A TETHER-FREE LIGHT SOLUTION

Experience the freedom of a wireless, tangle-free work environment with Atom Plus Wireless LED Light. This lightweight, versatile solution provides a bright, white spotlight that reaches 26,500 lux at max intensity. With magnetic clip batteries that attach easily and wirelessly without cables, you'll have a streamlined working environment, tangle-free. The Atom Plus wireless light is ideal for professionals who need a powerful lighting option for clinical or minor surgical procedures. The package comes with two rechargeable batteries, a charger, detachable filter, and screwdriver- all packed in a carrying case.

### VERSATILE AND ADAPTABLE

Recommended for clinical professionals performing procedures for:

- Facial Plastics
- Oral and Maxillofacial
- General Dentistry
- Cosmetic Dentistry

### MAGNETIC LOUPE ADAPTERS

With our selection of magnetic loupe adapters, quickly and easily fit your existing face shields or most major brands with the perfect accessory. The reliable magnet or attachment will provide all day hold.

### TOUCH ACTIVATION CONTROL

The Atom Plus batteries feature capacitive touch technology that can be used with gowns or gloves. The touch activation powers the light and controls three levels of brightness intensity providing you with even more flexibility during procedures or exams.

# TECHNICAL DATA

**Light Source:** LED refractive technology

**Lamp Illumination:** 26,500 lux

**Color Temperature:** 6,000k

**Illumination Spot Size:** 3.74 inches, fixed


**Construction:** Aluminum

**LED Light Module Weight:** 18 grams

**LED Warranty:** 5 years

**Battery:** Rechargeable Lithium Ion (charger included)

**Battery Warranty:** 1 year

<b>Battery Pack Options</b>	<b>Atom Plus Battery</b>	
<b>Battery Run-Time</b>	2 hours - 2 batteries included	
<b>Battery Pack Dimensions</b>	55mm x 28mm x 15mm	
<b>Battery Pack Weight</b>	23 grams	

## DESIGNED TO MEET YOUR NEEDS

Illuminate your workspace with Enova's incredible selection of innovative lighting solutions. Brighten deep surgical cavities with the Q-Series LED headlight or opt for the convenience of a cordless clip-on for minor procedures. Visit [www.goenova.com](http://www.goenova.com) or email customer service at [sales@goenova.com](mailto:sales@goenova.com) to learn more about custom lighting solutions for your practice.



**Qubit Cool LED Headlight**  
105,000 lux  
Tethered wearable battery  
8+ hour run-time



**Axis Light™**  
30,000 lux  
Wearable neck battery or traditional battery pack  
10+ hour run-time



**Atom Plus Wireless LED Light**  
26,500 lux  
Wireless Battery

To learn more about how Enova can help you see better, so you can perform better, visit [www.goenova.com](http://www.goenova.com), call 651.236.8857, or email [sales@goenova.com](mailto:sales@goenova.com).

### World Headquarters

Enova Illumination  
6305 Sandburg Road STE 100  
Minneapolis, MN 55427